


**NEW!!!**

# BD CARV II™ LX Laser/White Light Confocal

The best of both worlds



Helping all people  
live healthy lives

## Lamp and Lasers. No Limits.

Fill in your low light and high speed confocal applications with the new BD CARV II LX Confocal Imager. The BD CARV II LX combines both laser and white light illumination to satisfy a whole spectrum of wide field and confocal imaging applications.

**Discover more with BD – where confocal imaging technology knows no boundaries.**

**For more information visit us at [bdbiosciences.com/bioimaging](http://bdbiosciences.com/bioimaging)**

- Excitation light source combines Hg/metal halide lamp and up to 5 laser lines for ample power throughout the spectrum
- 3-5 fold increase in signal-to-noise at selected wavelengths
- Switch between wide field and confocal applications
- Suited for low-light applications
- Achieve fast confocal kinetic data